

Correction to "Protein Micropatterning in 2.5D

Citation for published version (APA):

van der Putten, C., Buskermolen, A. B. C., Werner, M., Brouwer, H. F. M., Bartels, P. A. A., Dankers, P. Y. W., Bouten, C. V. C., & Kurniawan, N. A. (2022). Correction to "Protein Micropatterning in 2.5D: An Approach to Investigate Cellular Responses in Multi-Cue Environments". *ACS Applied Materials & Interfaces*, 14(13), 15859. <https://doi.org/10.1021/acsami.2c04641>

Document license:
CC BY

DOI:
[10.1021/acsami.2c04641](https://doi.org/10.1021/acsami.2c04641)

Document status and date:
Published: 06/04/2022

Document Version:
Publisher's PDF, also known as Version of Record (includes final page, issue and volume numbers)

Please check the document version of this publication:

- A submitted manuscript is the version of the article upon submission and before peer-review. There can be important differences between the submitted version and the official published version of record. People interested in the research are advised to contact the author for the final version of the publication, or visit the DOI to the publisher's website.
- The final author version and the galley proof are versions of the publication after peer review.
- The final published version features the final layout of the paper including the volume, issue and page numbers.

[Link to publication](#)

General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal.

If the publication is distributed under the terms of Article 25fa of the Dutch Copyright Act, indicated by the "Taverne" license above, please follow below link for the End User Agreement:

www.tue.nl/taverne

Take down policy

If you believe that this document breaches copyright please contact us at:

openaccess@tue.nl

providing details and we will investigate your claim.

Correction to “Protein Micropatterning in 2.5D: An Approach to Investigate Cellular Responses in Multi-Cue Environments”

Cas van der Putten, Antonetta B. C. Buskermolen, Maike Werner, Hannah F. M. Brouwer, Paul A. A. Bartels, Patricia Y. W. Dankers, Carlijn V. C. Bouten, and Nicholas A. Kurniawan*

ACS Appl. Mater. Interfaces 2021, 13 (22), 25589–25598. DOI: 10.1021/acsami.1c01984



Cite This: *ACS Appl. Mater. Interfaces* 2022, 14, 15859–15859



Read Online

ACCESS |

Metrics & More

Article Recommendations

In the original version of this article, section *Materials and Methods*, Table 1 contains a mistake. In the description of “passivation with poly-L-lysine and mPEG-SVA”, the concentration of poly-L-lysine should be 0.1 mg/mL. In the revised table below, this correct concentration is included. This correction does not alter any conclusions of this work.

Table 1. Stepwise Description of Five Passivation Strategies for PDMS Surfaces

passivation strategy	description	duration and temp
gas-phase passivation with APTES	1. Place an open vial with APTES and the PDMS substrate in a hermetically sealed Petri dish.	1 h, room temp
	2. Gently wash three times with HEPES (0.1 M) (8 < pH < 8.5).	
	3. Incubate the PDMS substrate with 50 mg/mL mPEG-SVA in HEPES (0.1 M) (8 < pH < 8.5).	1 h, room temp
	4. Wash three times with PBS.	
liquid-phase passivation with APTES	1. Immerse the PDMS substrate in a 1% APTES solution in Milli-Q.	1 h, room temp
	2. Wash three times with PBS.	5 min/wash
	3. Wash three times with HEPES (0.1 M) (8 < pH < 8.5).	
	4. Incubate the PDMS substrate with 50 mg/mL mPEG-SVA in HEPES (0.1 M) (8 < pH < 8.5).	1 h, room temp
	5. Wash three times with PBS.	
vacuum passivation with APTES	1. Place the PDMS substrate and an open vial of APTES in a vacuum desiccator.	overnight, room temp
	2. Wash three times with HEPES (0.1 M) (8 < pH < 8.5).	
	3. Incubate the PDMS substrate with 50 mg/mL mPEG-SVA in HEPES (0.1 M) (8 < pH < 8.5).	1 h, room temp
	4. Wash three times with PBS.	
passivation with poly-L-lysine and mPEG-SVA	1. Pretreat the PDMS substrate with oxygen plasma at 20–100 W.	30–60 s
	2. Incubate the substrate with 0.1 mg/mL poly-L-lysine.	30 min, room temp
	3. Wash three times with HEPES (0.1M) (8 < pH < 8.5).	
	4. Incubate the PDMS substrate with 50 mg/mL mPEG-SVA in HEPES (0.1 M) (8 < pH < 8.5).	1 h, room temp
	5. Wash three times with PBS.	
passivation with Pluronic F-127	1. Immerse the PDMS substrate in a 1% Pluronic F-127 solution in PBS.	5 min
	2. Wash three times with PBS.	5 min/wash

Published: March 24, 2022

